

Mouse Bone Marrow Transplantation

A. Solutions

Flushing Solution (pH 7.3)

| | |
|---------|---|
| DMEM | 225 ml |
| FBS | 25 ml |
| heparin | 10 U/ml (3.125 ml of stock* per 250 ml) |

*heparin stock = 10 $\mu\text{g}/\mu\text{l}$
80 units heparin/mg = 0.08 units/ μg
So a 10 $\mu\text{g}/\mu\text{l}$ stock is 0.8 units/ μl .

ACK lysing buffer (pH 7.2-7.4)

| | |
|----------------------|---------|
| NH ₄ Cl | 4.15 g |
| KHCO ₃ | 0.5 g |
| Na ₂ EDTA | 18.6 mg |
| H ₂ O | 500 ml |

Washing Buffer

| | |
|---------|-------|
| 1X PBS | 50 ml |
| Albumin | 1g |

Suspending solution

| | |
|---------|--|
| DMEM | 50 ml |
| Albumin | 0.5 g |
| heparin | 5 U/ml (312.5 μl of stock*) |

B. Set up the day before the experiment

1. Prepare Flushing, Washing and Suspending solution fresh.
2. Sterilize dissecting tools.
3. Irradiate recipient mice.
4. Bring donor mice over from mouse house.

Take to mouse house:

University ID
P200 and sterile tips
mouse immobilizer
insulin syringes
250ml flask
cells
ice bucket