

### Protocol for Zymogram

#### Reagents:

1% Gelatin in H<sub>2</sub>O (Fisher Blood 275)

1% Casein (Sigma)

SDS-PAGE gel stock w/o urea

Wash buffer: 2.5% Triton X-100 in H<sub>2</sub>O (+0.02% NaN<sub>3</sub>)

Incubation Buffer: 50mM Tris-HCl (pH8.0), 5mM CaCl<sub>2</sub>, 0.02% NaN<sub>3</sub>

#### Gels:

Regular separating gel containing 10-12% substrate

Regular stacking gel

#### Protocol:

1. Collect media from cells (if desired, inactivate non-metalloproteases with PMSF and/or NEM)
2. Centrifuge to remove cellular debris (if necessary concentrate with centricon units or dialyze and lyophilize)
3. Add Laemmli loading buffer (**OMIT UREA AND REDUCING AGENTS, DO NOT HEAT**)
4. Load samples directly onto gel
5. Run gel
6. Wash 2X 20min in wash buffer
7. Wash 10min in incubation buffer
8. Place gel in sealable container with fresh incubation buffer and incubate at 37°C for 24h to 48h
9. Fix and stain with fresh Coomassie Blue solution
10. Destain with MeOH:AcOH:H<sub>2</sub>O(5:1:5)
11. Replace with 10% AcOH and continue destaining
12. Photograph and dry gel for storage