

Jax Generic Cre (Hotmaster)

Ingredient:	1 Reaction: (uL)	Master Mix Total (uL):	Sample Strain/ID Info:															
10X Buffer (Hotmaster)	3	0	<div style="border: 1px solid black; height: 100px; width: 100%;"></div> <table style="width: 100%; margin-top: 10px;"> <tr> <td style="width: 33%; text-align: center;">Sample</td> <td style="width: 33%; text-align: center;">Control</td> <td style="width: 33%; text-align: center;">Error</td> </tr> <tr> <td style="text-align: center;">Total:</td> <td style="text-align: center;">Total:</td> <td style="text-align: center;">Total:</td> </tr> <tr> <td style="background-color: yellow;"></td> <td style="background-color: yellow;"></td> <td style="background-color: yellow;"></td> </tr> <tr> <td colspan="3" style="text-align: center;">Total # of Sample Tubes:</td> <td style="border: 1px solid black; text-align: center;">0</td> </tr> </table>			Sample	Control	Error	Total:	Total:	Total:				Total # of Sample Tubes:			0
Sample	Control	Error																
Total:	Total:	Total:																
Total # of Sample Tubes:						0												
dNTPs (20mM)	0.3	0																
Cre-F (20uM)	1.5	0																
Cre-R (20uM)	1.5	0																
Int + Cont F (20uM)	1.5	0																
Int + Cont R (20uM)	1.5	0																
DMSO (2.5%)	0.5	0																
Taq Polymerase (Hotmaster)	0.2	0																
H ₂ O	18	0																

Add **28** uL master mix + **2** uL DNA (100ng/ul) per reaction tube

Total Volume: 30uL

Thermocycling:

94°C: 3 minutes

94°C: 30 seconds

51.7°C: 30 seconds

72°C: 1 minute

72°C: 2 minutes

10°C: ∞

x 35

Run on 3% gel for best resolution!

Expected Product Sizes:

Transgene: ~100bp

Internal Positive Control: 324bp

Primer Sequences:

IMR1084 (Cre-F) 5'-GCG GTC TGG CAG TAA AAA CTA TC-3'

IMR1085 (Cre-R) Rev 5'-GTC AAA CAG CAT TGC TGT CAC TT-3'

Int + Control F 5'- CTA GGC CAC AGA ATT GAA AGA TCT -3'

Int + Control R 5'-GTA GGT GGA AAT TCT AGC ATC ATC C-3'