

PRCE Genotyping Protocol

Minako K. Hashimoto-Partyka

(Hashimoto-Partyka MK, Lydon JP, and Iruela-ARispe ML. Genesis 2006)

For AB, AC, AD reactions:

<u>Reagent</u>	<u>Amount (for 30 μL reaction)</u>
H2O	12.3 μ L
10x Buffer (Hotmaster)	3.0 μ L
dNTP (2.5 mM)	2.5 μ L
Primer A	2.0 μ L
Primer B	2.0 μ L
Primer C	2.0 μ L
Primer D	2.0 μ L
Taq (Hotmaster)	0.2 μ L

For EC, ED reactions:

<u>Reagent</u>	<u>Amount (for 30 μL reaction)</u>
H2O	14.3 μ L
10x Buffer (Hotmaster)	3.0 μ L
dNTP 2.5 mM)	2.5 μ L
Primer E	2.0 μ L
Primer C	2.0 μ L
Primer D	2.0 μ L
Taq (Hotmaster)	0.2 μ L

Cycling conditions for both reactions:

1. 94°C for 2 min
2. 94°C for 20s
3. 59°C for 20s
4. 68°C for 40s
5. Repeat 2-4, 26x
6. 68°C for 5 min
7. 4°

Primer Sequences:

A – 5' TGT GCA CTT TTG GAG GCA AG 3'

B – 5' GTG GAG GCT TCT GGA CAG T 3'

C – 5' TAA AGC GCA TGC TCC AGA C 3'

D – 5' TGA TTT TGC CTT TGG CAG ATG 3'

E – 5' GGT CTC TGG CCT GAT TTT CC 3'

Primers

Genotype

AB PRCE, PR – ECKO(VEN), PRKO

AD PRKO