

Isolation of Ovarian Epithelial Cells

1. Remove ovaries from females
2. Label ___ eppendorfs
3. Add 300 μ l DMEM (w/o serum) to each
4. Add 3 μ l Dil-18* to all except the controls for each
5. Place an ovary in each tube
6. Short time point: Tubes _____
 - leave for ___ minutes in the 37°C incubator
 - fix the ovary from tube ___ and the control with 4%PFA for ___ hrs. (Eppendorf)
 - transfer ovary from tube ___ to 100 μ l DMEM
 - triturate
 - fill the rest of the tube with 4%PFA
7. Long time point: Tubes _____
 - leave for ___ hours in the 37°C incubator
 - fix the ovary from tube ___ and the control with 4%PFA for ___ hrs. (Eppendorf)
 - transfer ovary from tube ___ to 100 μ l DMEM
 - triturate
 - fill the rest of the tube with 4%PFA
8. After each time point:
 - Rinse triturated pieces and whole ovaries well with 1X PBS
 - Cut the whole ovaries in half on a flat surface
 - Return them to the PBS
 - Take one half to histology
 - Use the other for vibratome
9. Use a transfer pipette to transfer the triturated pieces to a slide
 - Remove excess PBS
 - If pieces are too big, cut with a razor blade
 - Add a couple of drops of TOPO3: Vectashield mounting media (1:500)
(Can also use 9 parts glycerol to 1 part PBS as mounting media)
 - View under confocal

*dioctadecyl-6,6'-di(4-sulfophenyl)-3,3,3',3'-tetramethylindocarbocyanine
(Molecular Probes)

Vibratome

1. 7% - 9% agarose (similar consistency to tissue)
 - Microwave – turn off as soon as it starts to boil
 - Set in ~38°C - 43°C water bath for 15 min.
2. Vibratome supplies:
 - Ice bucket
 - Chilled 1X PBS
 - Blade (cut in half)
 - Microslide (shelf above vibratome)
 - Cover slip forceps (flat head)
 - Vice adaptor
 - Baster
 - Kimwipes
 - Superglue (Can get at storeroom – IRU 200)
3. Lower the stage and tighten using knob facing you
4. Clean a Microslide
5. Label the molding block w/ sharpie to keep track of your ovaries
 - Keep the same orientation throughout
6. Transfer tissues from PBS to the microslide
7. Fill molding block w/ agarose ~3/4 full
8. Working quickly
 - place tissue in agarose
 - orient so the blade hits the narrower side first
 - leave on ice for ~1 minute to solidify
9. Slide out the solidified agarose
10. Trim excess agarose w/ blade
11. Fill the black tray w/ ice
12. Make sure blade angle is ~23-25°
 - avg. speed: 2-3
 - avg. amp.: 8-9
13. Superglue the small piece of agarose containing your tissue onto the metal tray
14. Place in front of fan to dry (only for ~5 min.)
15. Place tray in slot and tighten using knob on the right
16. Fill with chilled PBS
17. Advance (FFW) until you reach the tissue
 - slow down when you get close to the tissue
 - adjust speed/amp as needed
18. Keep the specimens separate in the tray as they're slivered off
19. Store in PBS or transfer directly to a microscope slide
 - use TOPO3:Mounting Media (1:500)

*dioctadecyl-6,6'-di(4-sulfophenyl)-3,3',3'-tetramethylindocarbocyanine
(Molecular Probes)