

## Media Recipes

A. Zovein 4-3-04

### MAEC media for Fetal liver:

Per 50ml media add:

Endothelial Mitogen (5mg/ml in PBS): 1ml

Heparin (salt): 50ul (1ul/ml)

VEGF 165 (100ng/ml): 25ul (50ng/ml)

MAEC above to make 2xabx

0.5ml P/S

50ul gent

### Pre-sort media:

25ml MAEC

25ml DMEM

0.75ml P/S

0.25ml noness aa

0.25ml l-glut

75ul gent

1.25ml HEPES buffer

### Post-sort media:

40ml MAEC

9ml FBS

0.5ml P/S

60ul gent

1.25ml HEPES

pre: 10%FBS, 2x abx, 10mmol HEPES

post: 40% FBS, 2x abx, 10mmol HEPES \*

\*10mmol HEPES: add 25ul per ml media

### Differentiation media:

Per 50ml media add:

50ul cAMP(25mg/ml)

15ul Retinoic acid (1mg/ml)

cAMP (protocol 0.5mM):

25ug/ml final []; stock make up to 25mg/ml

Retinoic Acid (protocol 1uM)

Stock 1mg/ml in DMSO

### 293T media

90ml DMEM

10ml FBS

1ml P/S