

## Bevan Lab standard internal solutions for electrophysiology

### K-MeSO<sub>4</sub> Internal Solution (chloride reversal = -83 mV)

	$\square$ (M)	MW	g/100 ml
<b>KCH<sub>3</sub>SO<sub>4</sub></b>	0.1300	150.2	1.9526
<b>NaCl</b>	0.0038	58.44	0.0222
<b>MgCl<sub>2</sub></b>	0.0010	203.31	0.0203
<b>HEPES</b>	0.0100	238.31	0.2383
<b>Phosphocreatine</b>	0.0050	453.4	0.2267
<b>Na<sub>4</sub>EGTA</b>	0.0001	468.3	0.0047
ATP and GTP are added on the day of the experiment			
<b>Na<sub>3</sub>GTP</b>	0.0004		20 $\mu$ l/ml of 20 mM stock
<b>Mg<sub>1.5</sub>ATP</b>	0.0020		20 $\mu$ l/ml of 100 mM stock

pH to 7.3 with 1 M KOH; dilute to 290 mmol/kg

For Ca<sup>2+</sup> imaging leave out Na<sub>4</sub>EGTA

### K-Gluconate Internal Solution-Bevan (chloride reversal = -83 mV)

	$\square$ (mM)	MW	g/100 ml
<b>K-Gluconate</b>	0.1400	234.25	3.2795
<b>NaCl</b>	0.0038	58.44	0.0222
<b>MgCl<sub>2</sub></b>	0.0010	203.31	0.0203
<b>HEPES</b>	0.0100	238.31	0.2383
<b>Na<sub>4</sub>-EGTA</b>	0.0001	468.3	0.0047
ATP and GTP are added on the day of the experiment			
<b>Na<sub>3</sub>GTP</b>	0.0004		20 $\mu$ l/ml of 20 mM stock
<b>Mg<sub>1.5</sub>ATP</b>	0.0020		20 $\mu$ l/ml of 100 mM stock

pH to 7.3 with 1 M KOH; dilute to 290 mmol/kg

**CsMeSO3 QX314**

	<b>□ (M)</b>	<b>MW</b>	<b>g/L</b>	<b>g/0.02913L</b>	<b>g/0.05826L</b>
<b>CsMeSO3</b>	0.1200	228	27.36	0.7970	1.5940
<b>NaCl</b>	0.0028	58.44	0.16	0.0048	0.0095
<b>HEPES</b>	0.0100	238.3	2.38	0.0694	0.1388
<b>TEA-Cl</b>	0.0050	165.7	0.83	0.0241	0.0483
<b>Na4EGTA</b>	0.0004	468.3	0.19	0.0055	0.0109
<b>QX314 HBr</b>	0.0050	343.31	1.72	0.0500	0.1000
<b>Phosphocreatine</b>	0.0050	453.4	2.27	0.0660	0.1321
<b>Spermine*</b>	0.0001	202.34	1 mL	29.12 µL	58.26 µL
ATP and GTP are added on the day of the experiment					
<b>Na3GTP</b>	0.0004			20 µl/ml of 20 mM stock	
<b>Mg1.5ATP</b>	0.0040			40 µl/ml of 100 mM stock	

pH to 7.3 with 1 M CsOH; dilute to 290 mmol/kg

**Spermine stock solution**

	<b>□ (M)</b>	<b>MW</b>	<b>g/L</b>	<b>g/100mL</b>	<b>g/10mL</b>
<b>* Spermine</b>	0.1	202.34	20.23	2.0234	0.2023

**CsCl QX314 Internal Solution**

	<b>□ (M)</b>	<b>MW</b>	<b>g/L</b>	<b>g/0.02913L</b>	<b>g/0.05826L</b>
<b>CsCl</b>	0.1350	168.4	22.73	0.6622	1.3245
<b>NaCl</b>	0.0036	58.44	0.21	0.0061	0.0123
<b>MgCl2.6H2O</b>	0.0010	203.31	0.20	0.0059	0.0118
<b>HEPES</b>	0.0100	238.3	2.38	0.0694	0.1388
<b>Na4EGTA</b>	0.0001	468.3	0.05	0.0014	0.0027
<b>QX314 HBr</b>	0.0100	343.31	3.43	0.1000	0.2000
ATP and GTP are added on the day of the experiment					
<b>Na3GTP</b>	0.0004			20 µl/ml of 20 mM stock	
<b>Mg1.5ATP</b>	0.0020			20 µl/ml of 100 mM stock	

pH to 7.3 with 1 M CsOH; dilute to 290 mmol/kg

**HEPES-buffered SIF**

	<b>□ (M)</b>	<b>MW</b>	<b>g/1000 ml</b>
<b>NaCl</b>	0.1400	58.44	8.1816
<b>Glucose</b>	0.0230	180.16	4.1437
<b>HEPES</b>	0.0150	238.31	3.5747
<b>KCl</b>	0.0030	74.55	0.2237
<b>MgCl<sub>2</sub></b>	0.0015	203.31	0.3050
<b>CaCl<sub>2</sub></b>	0.0016	147.02	0.2352

pH to 7.2 with 1 M NaOH; dilute to 300–310 mmol/kg