Immunoprecipitation Protocol for Tissue Samples

Reference

1. Paganoni, S., Bernstein, J., & Ferreira, A. (2010). Ror1-Ror2 complexes modulate synapse formation in hippocampal neurons. Neuroscience, 165(4), 1261–1274. https://doi.org/10.1016/j.neuroscience.2009.11.056

Solutions

- 1. RIPA Homogenization Buffer: 50 mM Tris-HCl pH 7.4, 150 mM NaCl, 1 mM EDTA
 - Add to 2 ml aliquot before use: 12 mM NaF, 2.7 mM Na₃VO₄, 3 mM PMSF, 1% Protease Inhibitor Cocktail
- 2. RIPA Lysis Buffer: 50 mM Tris-HCl pH 8, 150 mM NaCl, 0.1% SDS, 1% NP-40
 - Add to <u>2 ml aliquot</u> before use: 12 mM Sodium Deoxycholate (C₂₄H₃₉NaO₄), 119 mM NaF, 27 mM Na₃VO₄, 1% Protease Inhibitor Cocktail, 1% Triton X-100

Abbreviation

1. BSA: Bovine Serum Albumin

Protocol (All steps should be done at 4°C or on ice to minimize protease activity)

- 1. Homogenize tissue sample in 500 μ L of RIPA buffer in 1.5 mL Eppendorf tube.
- 2. Centrifuge @15,000 x g, 4°C for 30 min to pellet nuclei. The supernatant is the total cell lysate.
- 3. Transfer supernatant to fresh Eppendorf tube, add 15 μ L of a 1 mg/ml BSA stock solution (made in RIPA buffer) and rock for 1h at 4°C.
- 4. In fresh Eppendorf tube, combine supernatant with 10 μ L of primary antibody and incubate for 1-2 hours at 4°C rocking.
- 5. Add 25 μ L of pre-equilibrated Protein resin and incubate while rocking for 1 hour at 4°C.
 - a. Neutravidin comes in 50% v/v slurry, to pre-equilibrate in RIPA:
 - i. Aliquot ~26 µL slurry into Eppendorf tube
 - ii. Centrifuge at medium speed ~4,000 x g for 1 minute
 - iii. Decant supernatant
 - iv. Resuspend pellet in 13 μ L RIPA buffer and repeat washes ~3 times
 - v. After final spin, decant supernatant and resuspend pellet in 13 μ L RIPA
- 6. Spin at 14,000 x g for 1 minute to pellet, discard supernatant.
- 7. Wash at least 4 times with the lysing buffer (RIPA), centrifuge to collect pellet, discard the supernatant after each wash.
- 8. After final wash, discard the supernatant and resuspend pellet in 30-40 μ L of 2X Laemmli buffer.
- 9. Boil samples for 5-10 minutes to dissociate the complex.