

Membrane/Organelle Fractioning

References

1. Afreen, S., & Ferreira, A. (2019). Altered Cytoskeletal Composition and Delayed Neurite Elongation in tau₄₅₋₂₃₀-Expressing Hippocampal Neurons. *Neuroscience*, *412*, 1–15.
<https://doi.org/10.1016/j.neuroscience.2019.05.046>
2. Afreen, S., Riherd Methner, D. N., & Ferreira, A. (2017). Tau₄₅₋₂₃₀ association with the cytoskeleton and membrane-bound organelles: Functional implications in neurodegeneration. *Neuroscience*, *362*, 104–117.
<https://doi.org/10.1016/j.neuroscience.2017.08.026>
3. Nicholson, A. M., Wold, L. A., Walsh, D. M., & Ferreira, A. (2012). β -Amyloid carrying the Dutch mutation has diverse effects on calpain-mediated toxicity in hippocampal neurons. *Molecular medicine (Cambridge, Mass.)*, *18*(1), 178–185.
<https://doi.org/10.2119/molmed.2011.00366>

Solutions

1. Microtubule Stabilizing Buffer (MTSB): 130 mM HEPES, 4 mM MgCl₂, 10 mM EGTA, pH 6.9
2. Membrane/Organelle Fractioning Buffer: Microtubule Stabilizing Buffer + 0.02% Saponin

Abbreviations

1. MTSB: Microtubule Stabilizing Buffer
2. PBS: Phosphate Buffered Saline

Protocol

1. Before beginning, warm 10 ml of MTSB in 37°C in water bath.
2. Discard N2 media and wash once in PBS, remove and discard.
3. Add ~1 ml of detergent-free MTSB buffer for 30 seconds to stabilize, remove and discard.
4. Add ~1 ml of MTSB containing Saponin for 1 min to extract.
5. Rinse carefully in ~1 ml of detergent-free MTSB for 30 seconds, remove and discard.
6. Harvest fraction in 2X Laemmli buffer for Western blotting.
7. Boil samples for 10 minutes and store at -20°C.