

Sarkosyl Fractionation Assay

Cells or Recombinant Protein

References

1. Afreen, S., & Ferreira, A. (2022). The formation of small aggregates contributes to the neurotoxic effects of tau₄₅₋₂₃₀. *Neurochemistry international*, 152, 105252. <https://doi.org/10.1016/j.neuint.2021.105252>
2. Mocanu, M. M., Nissen, A., Eckermann, K., Khlistunova, I., Biernat, J., Drexler, D., Petrova, O., Schönig, K., Bujard, H., Mandelkow, E., Zhou, L., Rune, G., & Mandelkow, E. M. (2008). The potential for beta-structure in the repeat domain of tau protein determines aggregation, synaptic decay, neuronal loss, and coassembly with endogenous Tau in inducible mouse models of tauopathy. *The Journal of neuroscience : the official journal of the Society for Neuroscience*, 28(3), 737–748. <https://doi.org/10.1523/JNEUROSCI.2824-07.2008>
3. Chun, W., Waldo, G. S., & Johnson, G. V. (2007). Split GFP complementation assay: a novel approach to quantitatively measure aggregation of tau in situ: effects of GSK3beta activation and caspase 3 cleavage. *Journal of neurochemistry*, 103(6), 2529–2539. <https://doi.org/10.1111/j.1471-4159.2007.04941.x>
4. Cho, J. H., & Johnson, G. V. (2004). Glycogen synthase kinase 3 beta induces caspase-cleaved tau aggregation in situ. *The Journal of biological chemistry*, 279(52), 54716–54723. <https://doi.org/10.1074/jbc.M403364200>

Solutions

1. MES Buffer (store at -20°C): 20 mM MES, pH 6.8, 80 mM NaCl, 1 mM MgCl₂, 2 mM EGTA, 1% Protease Inhibitor Cocktail
2. Sarkosyl Buffer (store at -20°C): 500 mM NaCl, 10% (w/v) Sucrose, 1% (w/v) N-Lauroylsarcosine Sodium Salt (Sarkosyl)

Protocol

1. Day 1 - Aggregate protein following the normal procedure. This will yield protein at a final concentration of 4 μM. Set aside aliquots for 200 μL total that are aggregated simultaneously with those aggregated in the plate reader, but do not have ThioflavinS with them.
 - Note the presence of ThioflavinS at 100 μM has been shown to increase aggregation. We only use 20 μM ThioflavinS, but it is still best not to use the tau that has been used in the ThioflavinS aggregation assay.
2. Aggregate 5-6 hours at room temperature.
3. Add an equivalent volume of Ethanol to the unaggregated protein that has been kept on ice to the volume of Arachidonic Acid that was added to the aggregated protein.
4. Aliquot 20 μL of both the unaggregated and aggregated protein before spinning and mix 1:1 in 2X Laemmli Buffer. This will be total protein.
5. Spin all protein at 166,000 x g (30 psi) for 30 minutes.

6. Save supernatant (soluble protein). Bring final volume to 360 μ L with 2X Laemmli Buffer and boil for 10 minutes. Store at -20°C for WB analysis.
7. Re-suspend pellet (Insoluble protein) in 100 μ L of Sarkosyl Buffer.
 - Pellet is not visible so mark the “up” side of the tube.
 - This suspension should then be transferred to a 0.6 mL Eppendorf tube.
8. Vortex samples at room temperature for 30 minutes (tape to vortex).
9. Incubate overnight at 4°C with rocking.
10. Day 2 - Centrifuge the samples (Sarkosyl Soluble) for 1 hour at 166,000 x g.
11. Supernatant is soluble tau. Bring final volume to 360 μ L with 2X Laemmli and boil for 10 minutes.
12. Resuspend pellet (Sarkosyl Insoluble) in 360 μ L 2X Laemmli Buffer and boil for 10 minutes.
13. Analyze by Western blot.