

	Running Gels (30ml)				Stacking Gels (10ml)		
	8%	10%	12%	15%	3%	4%	5%
H2O	14	12	10	7	6.35	6.1	5.7
30% Acryl	8	10	12	15	1	1.33	1.66
4X Running Bf	7.5	7.5	7.5	7.5	--	--	--
Stacking Bf	--	--	--	--	2.5	2.5	2.5
10% SDS	0.3	0.3	0.3	0.3	0.1	0.1	0.1
10% APS	0.15	0.15	0.15	0.15	0.05	0.05	0.05
Temed	0.018	0.015	0.015	0.015	0.01	0.01	0.01

Wet Electrophoretic Transfer

Cut one sheet of nitrocellulose and four sheets of absorbent filter paper (Whatman 3MM) to the size of the gel.

Soak the nitro in distilled H₂O. Nitrocellulose should be wetted by carefully laying it on the surface of water. Allow the nitrocellulose to wet by capillary action (several minutes), then submerge it for 2 min. Move the membrane to soak in transfer buffer for 5 min. Wet the absorbent paper by soaking in transfer buffer.

Transfer Buffer 1: (proteins 20-400kD)

	conc.	for 1L
Tris	50mM	5.8g
Glycin	380mM	29g
SDS	0.10%	1g
MeOH	20%	200ml
Distilled H2O		to 1L

Transfer Buffer 2: (proteins <80kD)

	conc.	for 1L
	25mM	2.9g
	190mM	14.5g
	--	--
	20%	200ml
		to 1L