

## mRNA Synthesis

### Stocks:

5X Transcription buffer (supplied with enzyme)  
0.1M DTT (supplied with enzyme)  
ribonucleic acids 100mM each (Roche rNTP set)  
10mM rGTP dilution (from Roche rNTP set)  
5-methyl G (capG) 20mM (Roche 10904988001)  
RNasin (RNase inhibitor)  
RNA Polymerase  
Template at 0.5-1  $\mu\text{g}/\mu\text{l}$   
DEPC H<sub>2</sub>O

<b>Per Reaction:</b>	<b>25ul</b>	<b>50ul</b>
1X Transcription Buffer	5ul	10ul
10mM DTT	2.5ul	5ul
2mM each CTP, ATP, UTP	0.5ul (each)	1ul (each)
0.2mM GTP (1:10 stock)	0.5ul	1ul
1mM capG	1.25ul	2.5ul
1-2ug linearized template		
RNasin	1ul	2ul
Sp6 RNA Polymerase	1.25ul	2.5ul
DEPC H <sub>2</sub> O	to 25ul	to 50ul

Incubate at 37°C for 2hrs

DNase treat 30' @ 37°C                      1ul                      2ul

-run 1ul on 1% check gel (fast run on agarose gel, 5' or less)

-Bring volume to 100 $\mu\text{l}$  with DEPC H<sub>2</sub>O

-Extract 2X with 100 $\mu\text{l}$  Phenol-Chloroform

-Precipitate aqueous phase:

100ul RNA + 26ul 7.5M NH<sub>4</sub>OAc + 300ul 100% EtOH

-Spin 30' @ 4°C

-Remove supe and resuspend in 50ul DEPC H<sub>2</sub>O

-Precipitate again:

50ul RNA + 5ul 3M NaOAc + 200ul 100% EtOH

-Spin 30' @ 4°C

-Wash in 70% EtOH, spin 5' @ 4°C

-Resuspend in DEPC H<sub>2</sub>O. 30ul for a 25ul reaction is typical

-Injecting a 1:4 dilution (1ul in 4ul DEPC H<sub>2</sub>O) of this RNA (10nl) is about 500pg